
ViGO™ Energy Pal



user's manual

TABLE OF CONTENT

1
About Vigo Energy Pal

7
For Best Effect

2
Please Do Not

8
Positioning

3
What's in the Box

9
Detoxification Period

4
Switch On

10
Benefits to Expect

5
Bathed in Energies

11
Help Conditions

6
Stop Early

12
Warranty & Maintenance

13
Product Specifications

1

About Vigo Energy Pal

Vigo Energy Pal packs five key natural energies that mimic the harmonic powers from the earth itself in one device. The energies can activate body cells and enhance the human immune system, enabling the human body to operate at optimum levels and improving people's wellbeing.

“The body is self-healing, self-regulating, and self-regenerating.” Vigo Energy Pal essentially helps provide the energy needed to do this.

- Pulsed Electromagnetic Field (PEMF)
- High Voltage Electrostatic Potential (HELP)
- Far Infrared Ray
- Negative Ions
- Moxibustion

2

Please Do Not

Please do not use Vigo Energy Pal when you have

- a pacemaker implant
- an electro-cardio implant
- any electronic heart or lung apparatus

Please seek medical advice when you

- have infectious symptoms e.g. a fever
- are pregnant
- have a malignant tumor
- have a life-threatening illness

3

What's in the Box

What's in the box

- Vigo Energy Pal
- Carrying case with zip closure and leaflet storage flap
- Power cord
- Non-slip floor mat
- Testing ball in sealed plastic tube

start button

timer



internal fuse

power port

on/off switch



organizer

power code

non-slip floor mat

testing ball in tube



4

Switch On

Plug in the power cord securely both to Vigo Energy Pal and a power socket

Switch the machine on at the rear of the device (1).
The timer display (5) will show 60 minutes

Press the START button on the top of the device (2).
The lights will flash and the timer will commence a
countdown from 60

5

Bathed in Energies

Relax and enjoy your wellness therapy. Please be aware that it typically takes about 5 to 10 minutes for the surface temperature to reach its stable working value.

All functions and energies are automatically operated by the computerized programme for the duration of the 60-minute session.

When the timer has counted down from 60 minutes to zero, an alarm will sound for a few seconds.

6

Stop Early

You may stop the device at any time, by switching off at the rear of the device (1). For best results, it is recommended that Vigo Energy Pal is used for a minimum of 30 minutes daily.

The timer will reset the next time it is switched on, using the START button (2) on the top of the device will commence a countdown from 60.

7

For Best Effect

To build up a cumulative beneficial effect, it is suggested that Vigo Energy Pal is used for between 30 and 60 minutes every day. It is impossible to overdose, or use the device too much.

To get the very best from your device, it is recommended that you commence use by applying your feet or hands/palms onto the device treatment area (3), where the five natural energies are emitted.

After a few days, if you have a specific area of your body that requires more focus, move the device as near to the affected area as practical, for maximum therapeutic effect.

It is also recommended that you drink more water than usual, when using the device. This helps dispel toxins faster and facilitates blood and circulation improvements.

8

Positioning

The diagram on following website shows many various options. All are equally effective. It is just down to personal choice.

www.vigo-tech.com

9

Detoxification Period

The first couple of weeks will see a build up of cumulative beneficial effect in the body. However, certain new users have reported to go through a detoxification period, where some conditions can worsen slightly. This is completely normal and a necessary part of the healing process.

10

Benefits to Expect

- increased energy level
- improved sense of wellbeing
- reduced chronic pain
- less swelling
- better sleeping pattern
- accelerated injury recovery time
- sharper mental focus
- improved circulation
- enhanced immunity
- stabilised body functions e.g. blood pressure and blood sugar levels
- improved digestion processes
- less stressful and more relaxed

11

Help Conditions

There are countless case studies of people who have enjoyed astonishing results from using Vigo Energy Pal.

The most common reported improvements are:

- chronic pain
- circulatory problems
- hypertension
- gout digestion problems
- immunity difficulties
- asthma
- rhinitis
- skin complaints e.g. psoriasis
- insomnia
- oedema

12

Warranty & Maintenance

The only maintenance required, is simple cleaning of the surface area with a soft cloth, possibly infused with a non-toxic household disinfectant. Please do not use alcohol-based cleaners.

Warranty: 5 years

The Warranty will not cover damages caused by negligence, improper use, opening the chassis, or using the device for a purpose that it is not designed for.

13

Product Specifications

electrical	
power supply	220V/230V±10% , 50Hz±1Hz
power consumption	≤ 80W
fuse	F2AL, 1A / 250V, (F1AL, 2A / 250V) RP1-20, Ø5×2
timer	60 min
therapeutic	
pulsed electromagnetic field	
intensity	2000~ 20000 μTesla
frequency	50 Hz
working cycle	15 seconds
waveform	sin & cos constant
high electrostatic potential	
surface electrostatic potential	≥-100 KV/m
current	≤ 50 mA
frequency	Static/DC
far infrared	
peak wavelength range	9.0 ~ 9.9 μm
negative ion	
concentration	≥ 20000/cm ³
thermobustion	
temperature	37~43°C
physical	
dimensions (of the main unit)	390 mm × 350 mm × 95 mm
net weight	6 kg
chassis material	PC+ABS UL94_V0
load bearing capacity	max 135 kg
storage	
temperature	0°C ~ 55°C
humidity	≤ 80
atmospheric pressure	700 ~ 1060 hPa
maintenance	
installation	not required
maintenance	not required
warranty	
warranty	5 years

Read and fully understand this user's manual before
using Vigo Energy Pal.
Failure to follow operating instructions could result
in death or serious injury.

CONTACT

117 Siwei Road Yangmei
Taoyuan City, Taiwan 326

+886 (03) 425-9598
info@vigo-tech.com

www.vigo-tech.com